

Brucella Antibody - FOR EXPORT ONLY

SAMPLE REQUIRED:

Whole blood in a serum tube,
minimum of 2ml

BLOOD TUBE REQUIRED:

Serum tube, red top or gold top

Indications:

- For serological evaluation of dogs for export requirements only.
- Serology for suspected clinical cases of Brucella cannot be accepted by Vetnostics.

Collection Protocols:

1. A serum sample is collected into a red top or gold top tube.
2. The microchip number must be written on the sample tube along with the animal and owner names.
3. An AAHL request form must be included with the sample, along with our routine Vetnostics request form.